

Product Datasheet

Rat AFM ELISA Kit (orb437779)

Description	Rat Afamin (AFM) competitive inhibition ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P36953
Sensitivity	0.93ng/mL
Entrez	282708
Expiration Date	Please enquire.