

Product Datasheet

Canine aFP ELISA Kit (orb437781)

Description	Canine Alpha-Fetoprotein (aFP) Sandwich ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MJU5
Sensitivity	0.133ng/mL
Entrez	403551
Expiration Date	Please enquire.