

Product Datasheet

Equine aFP ELISA Kit (orb437782)

Description	Equine Alpha-Fetoprotein (aFP) Sandwich ELISA Kit
Reactivity	Equine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49066
Sensitivity	0.116ng/mL
Entrez	100034185
Expiration Date	Please enquire.