

Product Datasheet

Rat AGGF1 ELISA Kit (orb437797)

Description	Rat Angiogenic Factor With G Patch And FHA Domains 1 (AGGF1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.468-30ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F1MAK0
Sensitivity	0.19ng/mL
Entrez	310005
Expiration Date	Please enquire.