

Product Datasheet

Human AGRN ELISA Kit (orb437800)

Description	Human Agrin (AGRN) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	000468
Sensitivity	14.6pg/mL
Entrez	375790
Expiration Date	Please enquire.