

Product Datasheet

Human AGRP ELISA Kit (orb437801)

Description	Human Agouti Related Protein (AGRP) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	000253
Sensitivity	27pg/mL
Entrez	181
Expiration Date	Please enquire.