

Product Datasheet

Rat AGRP ELISA Kit (orb437803)

Description	Rat Agouti Related Protein (AGRP) Sandwich ELISA Kit
Reactivity	Rat
Range	9.375-600pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1MAG1
Sensitivity	3.5pg/mL
Entrez	25582
Expiration Date	Please enquire.