

Product Datasheet

Mouse AIF1 ELISA Kit (orb437818)

Description	Mouse AIF1 ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O70200
Sensitivity	12.2pg/mL
Entrez	11629
Expiration Date	Please enquire.