

Product Datasheet

Human AIMP1 ELISA Kit (orb437824)

Description	Human Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) Sandwich ELISA...
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q12904
Sensitivity	43pg/mL
Entrez	9255
Expiration Date	Please enquire.