

Product Datasheet

Mouse AIMP1 ELISA Kit (orb437825)

Description	Mouse Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) Sandwich ELISA...
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q3UZG4
Sensitivity	14.8pg/mL
Entrez	13722
Expiration Date	Please enquire.