

Product Datasheet

Mouse AIRE ELISA Kit (orb437827)

Description	Mouse Autoimmune Regulator (AIRE) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z0E3
Sensitivity	0.058ng/mL
Entrez	11634
Expiration Date	Please enquire.