

Product Datasheet

Human AK3 ELISA Kit (orb437828)

Description	Human Adenylate Kinase 3 (AK3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9UIJ7
Sensitivity	0.146ng/mL
Entrez	50808
Expiration Date	Please enquire.