

Product Datasheet

Bovine aLA ELISA Kit (orb437829)

Description	Bovine Alpha-Lactalbumin (aLA) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	milk and other biological fluids.
Sample Types	milk and other biological fluids.
Uniprot ID	P00711
Sensitivity	0.32ng/mL
Entrez	281894
Expiration Date	Please enquire.