

Product Datasheet

Human ALAD ELISA Kit (orb437830)

Description	Human Aminolevulinate Delta Dehydratase (ALAD) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, erythrocyte lysates and other biological fluids.
Sample Types	serum, erythrocyte lysates and other biological fluids.
Uniprot ID	P13716
Sensitivity	0.128ng/mL
Entrez	210
Expiration Date	Please enquire.