

Product Datasheet

Rat ALAD ELISA Kit (orb437831)

Description	Rat Aminolevulinate Delta Dehydratase (ALAD) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P06214
Sensitivity	0.128ng/mL
Entrez	25374
Expiration Date	Please enquire.