

Product Datasheet

Human aLA ELISA Kit (orb437832)

Description	Human Alpha-Lactalbumin (aLA) Sandwich ELISA Kit
Reactivity	Human
Range	6.25-400pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	breast milk and other biological fluids.
Sample Types	breast milk and other biological fluids.
Uniprot ID	P00709
Sensitivity	2.58pg/mL
Entrez	3906
Expiration Date	Please enquire.