

Product Datasheet

Rat ALAS1 ELISA Kit (orb437833)

Description	Rat Aminolevulinate Delta Synthase 1 (ALAS1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P13195
Sensitivity	0.36ng/mL
Entrez	65155
Expiration Date	Please enquire.