

Product Datasheet

Canine ALB ELISA Kit (orb437834)

Description	Canine Albumin (ALB) competitive inhibition ELISA Kit
Reactivity	Canine
Range	1.56-100µg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P49822
Sensitivity	0.49µg/mL
Entrez	403550
Expiration Date	Please enquire.