

Product Datasheet

Human ALDH1A1 ELISA Kit (orb437839)

Description	Human Aldehyde Dehydrogenase 1 Family, Member A1 (ALDH1A1) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00352
Sensitivity	0.236U/L
Entrez	216
Expiration Date	Please enquire.