

Product Datasheet

Human ALDH1A2 ELISA Kit (orb437841)

Description	Human Aldehyde Dehydrogenase 1 Family, Member A2 (ALDH1A2) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	094788
Sensitivity	0.197ng/mL
Entrez	8854
Expiration Date	Please enquire.