

Product Datasheet

Rat ALDOA ELISA Kit (orb437848)

Description	Rat Aldolase A, Fructose Bisphosphate (ALDOA) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05065
Sensitivity	0.058ng/mL
Entrez	24189
Expiration Date	Please enquire.