

Product Datasheet

Rat ALDOB ELISA Kit (orb437850)

Description	Rat Aldolase B, Fructose Bisphosphate (ALDOB) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q66HT1
Sensitivity	37pg/mL
Entrez	24190
Expiration Date	Please enquire.