

Product Datasheet

Human Alk-Smase ELISA Kit (orb437855)

Description	Human Alkaline Sphingomyelinase (Alk-Smase) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6UWV6
Sensitivity	1.39ng/mL
Entrez	339221
Expiration Date	Please enquire.