

Product Datasheet

Mouse Alk-Smase ELISA Kit (orb437856)

Description	Mouse Alkaline Sphingomyelinase (Alk-Smase) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q3TIW9
Sensitivity	0.051ng/mL
Entrez	238011
Expiration Date	Please enquire.