

Product Datasheet

Mouse ALPL ELISA Kit (orb437869)

Description	Mouse Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P09242
Sensitivity	0.057ng/mL
Entrez	11647
Expiration Date	Please enquire.