

Product Datasheet

Rat ALPL ELISA Kit (orb437870)

Description	Rat Alkaline Phosphatase, Liver/Bone/Kidney (ALPL) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08289
Sensitivity	0.31ng/mL
Entrez	25586
Expiration Date	Please enquire.