

Product Datasheet

Human ALPP ELISA Kit (orb437872)

Description	Human Alkaline Phosphatase, Placental (ALPP) Sandwich ELISA Kit
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05187
Sensitivity	5.1pg/mL
Entrez	250
Expiration Date	Please enquire.