

Product Datasheet

Human aMACR ELISA Kit (orb437882)

Description	Human Alpha-Methylacyl Coenzyme A Racemase (aMACR) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UHK6
Sensitivity	0.059ng/mL
Entrez	23600
Expiration Date	Please enquire.