


Product Datasheet

DDDDK tag antibody (orb43789)

Description	Mouse Monoclonal to DDDDK tag.
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide: DYKDDDDK conjugated to KLH
Target	DDDDK tag
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Expiration Date	12 months from date of receipt.



Immunocytochemistry
(confocal
microscopy...)