

Product Datasheet

Human AMY1 ELISA Kit (orb437897)

Description	Human Amylase Alpha 1, Salivary (AMY1) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva and other biological fluids.
Sample Types	serum, plasma, saliva and other biological fluids.
Uniprot ID	P04745
Sensitivity	1.32ng/mL
Entrez	276
Expiration Date	Please enquire.