

Product Datasheet

Bovine ANGPT1 ELISA Kit (orb437904)

Description	Bovine Angiopoietin 1 (ANGPT1) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O18920
Sensitivity	0.195ng/mL
Entrez	282140
Expiration Date	Please enquire.