

Product Datasheet

Canine ANGPT1 ELISA Kit (orb437905)

Description	Canine Angiopoietin 1 (ANGPT1) Sandwich ELISA Kit
Reactivity	Canine
Range	2.343-150ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q60FC1
Sensitivity	0.93ng/mL
Entrez	403656
Expiration Date	Please enquire.