

Product Datasheet

Rat ANGPTL3 ELISA Kit (orb437922)

Description	Rat Angiopoietin Like Protein 3 (ANGPTL3) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q5I0L8
Sensitivity	6.4pg/mL
Entrez	502970
Expiration Date	Please enquire.