

Product Datasheet

Rat ANKRD1 ELISA Kit (orb437933)

Description	Rat Ankyrin Repeat Domain Protein 1 (ANKRD1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8R560
Sensitivity	0.058ng/mL
Entrez	27064
Expiration Date	Please enquire.