

Product Datasheet

Canine ANP ELISA Kit (orb437937)

Description	Canine Atrial Natriuretic Peptide (ANP) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07499
Sensitivity	6.0pg/mL
Entrez	608289
Expiration Date	Please enquire.