

Product Datasheet

Mouse ANXA2 ELISA Kit (orb437949)

Description	Mouse Annexin A2 (ANXA2) Sandwich ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07356
Sensitivity	0.61ng/mL
Entrez	12306
Expiration Date	Please enquire.