

Product Datasheet

Human ANXA4 ELISA Kit (orb437952)

Description	Human Annexin A4 (ANXA4) Sandwich ELISA Kit
Reactivity	Human
Range	18.75-1200pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P09525
Sensitivity	6.9pg/mL
Entrez	307
Expiration Date	Please enquire.