

Product Datasheet

Mouse ANXA5 ELISA Kit (orb437955)

Description	Mouse Annexin A5 (ANXA5) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P48036
Sensitivity	5.8pg/mL
Entrez	11747
Expiration Date	Please enquire.