

Product Datasheet

Monkey ANXA5 ELISA Kit (orb437957)

Description	Monkey Annexin A5 (ANXA5) Sandwich ELISA Kit
Reactivity	Monkey
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q4R4H7
Sensitivity	0.128ng/mL
Entrez	101926167
Expiration Date	Please enquire.