

Product Datasheet

Human AOA^H ELISA Kit (orb437960)

Description	Human Acyloxyacyl Hydrolase (AOAH) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28039
Sensitivity	0.059ng/mL
Entrez	313
Expiration Date	Please enquire.