

Product Datasheet

Bovine AOX1 ELISA Kit (orb437962)

Description	Bovine Aldehyde Oxidase 1 (AOX1) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P48034
Sensitivity	0.107ng/mL
Entrez	338074
Expiration Date	Please enquire.