

Product Datasheet

Human AP13 ELISA Kit (orb437964)

Description	Human Apelin 13 (AP13) competitive inhibition ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	37.2pg/mL
Expiration Date	Please enquire.