

Product Datasheet

Rat APAF1 ELISA Kit (orb437972)

Description	Rat Apoptotic Peptidase Activating Factor 1 (APAF1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9EPV5
Sensitivity	0.062ng/mL
Entrez	78963
Expiration Date	Please enquire.