

Product Datasheet

Mouse APLN ELISA Kit (orb437982)

Description	Mouse Apelin (APLN) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9R0R4
Sensitivity	37.5pg/mL
Entrez	30878
Expiration Date	Please enquire.