


Product Datasheet

GFAP antibody [Out of stock] (orb43799)

Description

Mouse Monoclonal to GFAP.


Reactivity


Feline, Human, Porcine

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications

Immunogen

Pellet of porcine brain cold-stable proteins after depolymerization of microtubules.

Immunohistochemistry staining of human c...

Target

GFAP

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration


1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only



Immunohistochemistry staining of human b...

Application notes

Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: human brain (cortex, cerebellum). The antibody GF-01 strongly stains astrocytes in human brain tissue sections but it is essentially negative on mouse and rat tissues. Immunocytochemistry: Recommended dilution: 5-10 µg/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P14136

Entrez

2670

Dilution Range

Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: human brain (cortex, cerebellum). The antibody GF-01 strongly stains astrocytes in human brain tissue sections but it is essentially negative on mouse and rat tissues. Immunocytochemistry: Recommended dilution: 5-10 µg/ml.

Expiration Date

12 months from date of receipt.