

Product Datasheet

Monkey APOA1 ELISA Kit (orb437994)

Description	Monkey Apolipoprotein A1 (APOA1) Sandwich ELISA Kit
Reactivity	Monkey
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	F7FQM7
Sensitivity	0.62ng/mL
Entrez	696969
Expiration Date	Please enquire.