

Product Datasheet

Human APOA2 ELISA Kit (orb437995)

Description	Human Apolipoprotein A2 (APOA2) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02652
Sensitivity	25.4ng/mL
Entrez	336
Expiration Date	Please enquire.