

Product Datasheet

Mouse APOA2 ELISA Kit (orb437996)

Description	Mouse Apolipoprotein A2 (APOA2) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09813
Sensitivity	28ng/mL
Entrez	11807
Expiration Date	Please enquire.