

Product Datasheet

Porcine APOA2 ELISA Kit (orb437997)

Description	Porcine Apolipoprotein A2 (APOA2) Sandwich ELISA Kit
Reactivity	Porcine
Range	4.687-300µg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1S1A9
Sensitivity	1.2µg/mL
Entrez	100153243
Expiration Date	Please enquire.