

Product Datasheet

Human APOA5 ELISA Kit (orb438003)

Description	Human Apolipoprotein A5 (APOA5) Sandwich ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6Q788
Sensitivity	26.73pg/mL
Entrez	116519
Expiration Date	Please enquire.